

Welcome to Chem 101 Lab

Introduction & Stoichiometry Experiment



TA Contact Info

❑ Name: Ellie Zheng

❑ Department of Chemistry, Duke University

❑ Email: lianjun.zheng@duke.edu

Lab Syllabus Highlights

- ❑ Lab is worth **25%** of your grade
- ❑ Type of Assignments
 - Lab Report (Lab Notebook pages) due **1 week** after completing the experiment
 - Due to the **Undergraduate Drop-box, Room 1219:**
- ❑ Penalty for Late assignments – 10% first day then 5% per day
- ❑ Grades will be posted on Sakai on the Lab Site

Lab Syllabus Highlights

☐ Attendance Policy

- **NO** lab make-ups
- **ONLY** one Excused absence, submit the Excused Absence for Short Term Illness form (to Lab manager and me)


☐ **Individual** reporting

- Pair work: same data, different discussion

☐ Tardiness = 1 point deducted from preparedness grade

☐ Arrive after presentation = sent home, **0 point** for that day's lab

Graded Assignments



Notebook Pages (30 points)

All of the following sections should be included:

(Fill out BEFORE lab):

1. **Name, course and section number, experiment title, and date**
2. **Objective** – Statement explaining the purpose of the experiment
3. **Procedure** – outline of the steps you will perform, organized and concise

(Fill out DURING lab):

4. **Observations** – Record your observations & measurements, procedural changes, etc.
5. **Calculations** – Show all work for your calculations and include any graphs
6. **Results** – Record in table format; these are the final results of the experiment
7. **Discussion** – Complete after Lab; Write in paragraph form. Demonstrate your understanding of the experiment, summarize your results and state your confidence in them, and discuss sources of error in the experiment, if applicable.

*More details could be found on Sakai

*Always follow the instructions

Notebook Pages

- left– Procedure
- right–Observation & Measurement

Or

- Leaving blank between procedures for observation;
 - Put measurement tables/record in the end
- Depends on how much observation/measurement to record

CONFIDENTIAL

109

Outline of Procedure

1. _____
2. _____
3. _____
4. _____

Observation & Measurement

- _____
- _____

Table1 Mass of XX

Or

Step 1: %&^*^(((&&%^%\$%

Corresponding observation

Step 2: (&*##&(_^\$*(

Corresponding observation

Measurements:

CONFIDENTIAL

Graded Assignments



Sapling Pre-lab Assignments – (15 points) Found on Sapling Learning; Due by 8:30 or 1:25 pm the day of your lab meeting.

Balance Czar – (10 points) Given for cleaning the balances and balance area for one experiment – **once a semester**

Cleanliness Czar – (10 points) Given for cleaning common areas, cleaning and returning shared equipment to proper place – **once a semester**

Bin Check – (10 points) Your hood randomly checked by TA after lab to make sure each piece of glassware is clean and in the correct bin – **once a semester**

Preparedness & Technique Grade – (10 points) Grade assigned as an evaluation of your overall performance in lab

Balance & Cleanliness Czar

Hood#	Name	Email	Balance Czar	Cleanliness Czar
1			Sept. 20	Oct. 4
			Sept. 20	Oct. 4
2			Oct. 4	Sept. 20
			Oct. 4	Sept. 20
3			Oct. 18	Nov. 8
			Oct. 18	Nov. 8
4			Nov. 8	Nov. 22
			Nov. 8	Nov. 22
5			Nov. 22	Oct. 18
			Nov. 22	Oct. 18

Writing Sessions

- Designed to give you a basic introduction to scientific writing and the various types of scientific literature.
- Skills you can apply to any science class not only Chemistry.
- You will meet for three 1-hr sessions on the weeks listed in your schedule in Room 1243.
- Attendance/Participation is required.
- First Session is next week – Watch the Announcements on Sakai for more information.

Laboratory Safety

In the case of an Accident (Broken Glass, Fire, Chemical Spill):

- ✓ **ALWAYS** Notify the TA
- ✓ Know where the following Safety Equipment is located:




Laboratory Safety

In the case of an Accident (Broken Glass, Fire, Chemical Spill):

- ✓ Dispose of waste, broken glass and syringes in the appropriate places!



- ✓ Know the Fire Evacuation Route



- ✓ Without permission, Never Pour Chemicals Down the Drain

Laboratory Safety

✓ Always wear proper lab ATTIRE!



ALWAYS wear
Safety Glasses



Wear Long pants



NO open toed shoes



ALWAYS wear gloves



Lab Partner/Fume Hood

- ☐ After Pre-Lab discussion, pick a lab partner.
- ☐ Each pair of students will be assigned to 1 fume hood. This is your hood for the rest of the semester
- ☐ In the cabinets under each hood, there are 4 glassware bins
- ☐ Students are responsible for checking the contents of each bin and ensure that the glassware is CLEAN and DRY after each experiment!
- ☐ If something is missing or broken, it can be replaced from the stockroom (**Room 1237**). Extra equipment can be placed in the bin at the front of the room.

Today's Experiment - Stoichiometry

Experimental Objective/Challenge – Apply several fundamental chemical concepts and work as a team to collect and analyze data/observations about a chemical reaction.

Background Concepts:

- 1) Conservation of Mass
- 2) Law of Definite Proportions
- 3) Dalton's Atomic Theory

Reaction - $A + B \longrightarrow ?$

What we know about the reaction:

- Reagents A & B combine to form a white, fluffy solid – a precipitate
- Without knowing the reagents identity, we still know:

Mass of the precipitate = the amount of product formed



Two Different Trials

Trial A – How much product forms when you change the concentration of reagent A?

- Add 10 mL of 0.25, 0.5 M, 0.75 M, and 1 M Reagent A to 10 mL of 1 M Reagent B

Trial B – How much product forms when you vary the volume of reagent A added to the reaction?

- Add 8, 10, 12, 14 mL of 1 M Reagent A to 10 mL of Reagent B



Hood# 1~2: Trial A

Hood# 3~5: Trial B

Then we'll share the data within the class.

Questions to Consider:


- Predict what your graph will look like.
(Volume/Concentration of A vs. Mass of your product)




- How to determine which reagent is in excess?

Experiment Details



- We will use **Vacuum Filtration** to separate the product and then weigh it to find its mass using the Balance



Experiment Details



Graduated Cylinder



**Adding drops
& Rinsing**



Spatula



Plastic Pipet

Experiment Details



Watch glass



Evaporating dish



Experiment Details

- **Transferring**

- From bottle to beaker: pour directly, **roughly**
- From beaker to graduate cylinder: pour at first (like, until 9 mL), then using pipet to drop till 10.0 mL, **precisely**

** Using the proper size of graduate cylinder

- **Rinsing**

- Using pipet or squirt bottle



Experiment Details



- **Take masses** of filter paper and watch glasses, **label** beakers and graduated cylinders before use
- Set up **TWO** vacuum filtration at the same time. Turn on the vacuum just before transferring.
- Think about probable sources of error:
- **Be sure to:**
 - Rinse solid adhering to your spatula back into test tubes
 - Rinse the test tube to transfer your solid completely
 - Let each reaction mixture sit for full ~30 min (?)
 - Break up the precipitate to completely dry
 - Dry over vacuum for 15 min

loss

add

Safety & Waste

- Reagent B is an eye and skin irritant; Always wear safety glasses and rinse skin immediately if there is contact!
- Solids and filter paper should be disposed of in the Solid Waste Container.
- All excess reagents should be placed in Waste Containers.
- ONLY the reaction solution can go down the drain

Before Leaving the Lab

- Once you have all of your data collected, enter it into my excel spreadsheet.
**For Trial B, check the volume of reagent A before entering it.
- After I check the cleanliness of your hood and sign your notebook, you're free to go.
** Hood#1 ~ Balance czar, Hood#2 ~ Cleanliness czar

Note:

- (1)Data spreadsheet, (2) Detailed instruction for Reporting & Discussion, (3) Rosters and (4) Balance & Cleanliness Czar will be sent to you via email later today.
- You'll have to discuss both trial A & B. Using the whole data we collected in the spreadsheet.